

Knotless/Crimpleless Lateral Stabilization System

Ti Posts

OrthoZip

Isometric. Adjustable. Strong.

From the engineers that delivered the first crimping system, snap-off bone anchors, cruciate needles, swaged-on monofilament and custom instruments.




Features & Benefits:

- Strong, Fast & Effective
- Placed at Isometric Sites
- More Uniform Tension
- Smooth Rotation Posts
- Reduced Surgery Time
- **Less Inventory**
- Adjustable to Desired Size
- **Out Performs** Monofilament and Other Fibers*
- Many Applications
- Versatile
- Cost Effective
- Mechanically Tested
- High-Strength
- Minimal Elongation
- Available in USP #2 or #5
- Possible Minimally Invasive surgery

*See data on reverse side

OrthoZIP

KNOTLESS/CRIMPLESS LATERAL STABILIZATION SYSTEM



ORTHOZIP PRODUCTS

OrthoZip Rotax Posts


- Titanium posts provide smooth, stress-free surface for suture to rotate on
- Post heads keeps suture securely on post body
- Self tapping, Hex drive
- Available in many lengths

OrthoZip

- Made from UHMWPE (Ultra High Molecular Weight Polyethylene) orthopedic suture
- Available in USP #2 and #5
- >2X the strength of standard fiber loops
- Blue tracer fiber for easy visibility
- Pre-sterile ready to use
- Many applications

Lateral Stabilization

- High grade medical certified materials
- Very high strength with minimal elongation
- No knot, No crimp, No expensive instrumentation



*Fiberwire is a registered trademark of Arthrex. Orthofiber is a registered trademark of MWL. **Data collection for monofilament, Fiberwire & Orthofiber is from "Mechanical testing of orthopedic suture materials used for extra-articular stabilization of canine cruciate ligament-deficient stifles College of Veterinary Medicine, University of Minnesota."

EVEROST #	PATTERSON #	DESCRIPTION	QTY
005630	07-891-8486	OrthoZip Adjustable Locking Loop, UHMWPE # 2	
005632	07-891-8487	OrthoZip Adjustable Locking Loop, UHMWPE # 5	
005631	07-891-8876	OrthoZip Adjustable Locking Loop, UHMWPE # 2, Box of 6	
005633	07-891-8877	OrthoZip Adjustable Locking Loop, UHMWPE # 5, Box of 6	
263102	07-891-8873	Kit, OrthoZip Starter 3.5mm	
		OrthoZip Adjustable Locking Loop, UHMWPE # 2, Box of 6	
		3.5mm x 18mm Rotax OrthoZip Post	1
		3.5mm x22mm Rotax OrthoZip Post	2
		3.5mm x 24mm Rotax OrthoZip Post	2
263101	07-891-8874	Kit, OrthoZip Starter 2.7mm	
		OrthoZip Adjustable Locking Loop, UHMWPE # 2, Box of 6	1
		2.7mm x 14mm Rotax OrthoZip Post	2
		2.7mm x 18mm Rotax OrthoZip Post	2
		2.7mm x22mm Rotax OrthoZip Post	2
263100	07-891-8875	Kit, OrthoZip Starter 2.4mm	
		OrthoZip Adjustable Locking Loop, UHMWPE # 2, Box of 6	1
		2.4mm x 12mm Rotax OrthoZip Post	2
		2.4mm x 16mm Rotax OrthoZip Post	2
		2.4mm x20mm Rotax OrthoZip Post	2



800-225-7911 | www.pattersonvet.com | EXCLUSIVELY DISTRIBUTED BY



855.EVEROST(383 7678) • info@everost.com • www.everost.com